



Komet Medical[®] pins for knee joint endoprosthesis.















Once again, our new product range for knee joint endoprosthesis has shown that Komet Medical is a full-range supplier of motor-driven tools for surgical power systems. Pins play a decisive role in the fixation of knee templates. The Komet Medical range comprises different types of pins, all of which are available from stock in sterile or non-sterile packaging.



- Square
- Hexagon
- Roundshank



Range | Pins

	GP003.070.023	(Ø = mm)	L = mm	KL/HL= mm	Schaft	Wrench width (mm)
Pin		2,3	70	-	Roundshank	-
	GP002.089.032	(Ø = mm)	L = mm	KL/HL= mm	Shank	Wrench width (mm)
Pin		3,2	89,5	83,0	Square	2,15
	GP001.089.032	(Ø = mm)	L = mm	KL/HL= mm	Shank	Wrench width (mm)
Universal Pin		3,2	89	70,0	Square	2,25
	GS905.030.032	(Ø = mm)	L = mm	KL/HL= mm	Shank	Wrench width (mm)
Screw		3,2	30	26,0	Hexagon	3,5
	GS901.033.032	(Ø = mm)	L = mm	KL/HL= mm	Shank	Wrench width (mm)
Screw		3,2	33	30,0	Hexagon	2,5
	GS901.050.032	(Ø = mm)	L = mm	KL/HL= mm	Shank	Wrench width (mm)
Screw		3,2	50	46,0	Hexagon	2,5
	GS905.050.032	(Ø = mm)	L = mm	KL/HL= mm	Shank	Wrench width (mm)
Screw		3,2	50	46,0	Hexagon	3,5
	GS903.062.032	(Ø = mm)	L = mm	KL/HL= mm	Shank	Wrench width (mm)
Screw		3,2	62	46,0	Square	2,15
	GS902.150.032	(Ø = mm)	L = mm	KL/HL= mm	Shank	Wrench width (mm)
Long Pin		3,2	150	51,0	Roundshank	-
	GS904.100.030	(Ø = mm)	L = mm	KL/HL= mm	Shank	Wrench width (mm)
Long Pin		3,0	100	26,0	Square	2,15
	GS904.150.030	(Ø = mm)	L = mm	KL/HL= mm	Shank	Wrench width (mm)
Long Pin		3,0	150	26,0	Square	2,15
	GS904.100.040	(Ø = mm)	L = mm	KL/HL= mm	Shank	Wrench width (mm)
Long Pin		4,0	100	25,0	Square	2,15
	GS900.100.000				L = mm	Shank
Pin Adapter (for hexagonal screws and pins up to wrench width 2,5)					100	Zimmer/Hudson
	GS900.130.000				L = mm	Shank
Pin Adapter (for hexagonal screws and pins up to wrench width 2,5)					130	Zimmer/Hudson



Square



Hexagon



Roundshank